
Our Product Introduction

### Long Pulse Nd YAG Laser Hair Removal Machine Portable For Commercial

#### **Basic Information**

. Place of Origin: Shandong, China Brand Name: Healthy&Beauty Model Number: Picosecond Laser \$1,380.00/sets 1-99 sets • Packaging Details: outside:carton box



#### **Product Specification**

Name: Long Pulse Nd YAG Laser Hair Removal

inside:aluminium box

Machine

NO · Q-Switch:

. Laser Type: Nd: Yag Laser Style: Portable Type: Laser

Acne Treatment, Anti-Puffiness, Pigment Feature:

Removal, Pigmentation Correctors, Tattoo

Removal

· Application: For Commercial

 After-sales Service Free Spare Parts, Online Support, Video Provided: Technical Support, Field Installation,

Commissioning And Training, Field Maintenance And Repair Service

· Warranty: 2 Years

· Colour: Red,blue,slivery,pink,white,and So On


OEM ODM Services:

# TREATMENT EFFECT \_tattoo Evebrow Tattoo Age Pigment

#### More Images









#### **Product Description**

Air Water Temperature Control Nd Yag Laser Tattoo Removal Machine For Big Promotion

#### **Products Description**

## Picosecond Facial Rejuvenation active Q-Switched Tattoo Removal nd yag laser




#### Treatment Theropy:


Principle of laser tattoo removal: 1. Laser tattoo removal is the only method proven to remove tattoos without scarring. 2. To remove a tattoo, the laser is pulsed over the tattoo, directing its light energy towards the ink. 3. The energy is absorbed by the ink particles which then shatter into tiny fragments. 4. In the days following the laser procedure, the body's immune system will flush away the shat- tered ink particles and cause the tattoo to fade. 5. With more treatments, more ink will be shattered leaving the skin free of the tattoo.


#### Applications for picosecond laser:


1. Tattoo removal 2. Eyebrow embroidery removal 3. Eyeliner removal 4. Remove chloasma 5. Skin whitening 6. Removal of birthmark 7. Onychomycosis 8. Excision of nevus 9. Senile plaque removal 10. Melasma removal

### Principle


















#### Portable picosecond laser advantages High energy & Shot pulse duration

The power emit in trillions of a second to skin tissue, which can explore the pigment into much more small pieces than traditional Q:switch laser. The pigment pieces can be absorbed more easily and faster.

**Minimal risk and less pain** Since the pulse duration has been reduced thousand of time, less damage of skin is realized. The treatment is more safe and less pain.

**No anesthesia during the treatment, No inflammation** Fast For tattoo removal, usually 6-8 treatment is necessary by Traditional Q:switch laser. With pico second laser, only 3-4 treatment is needed only for completely removal.

**High performance** Completely removal of tattoo and pigment can be reached with Pico second technology. Suitable for skin color from I-IV skin Type,no hype pigmentation after treatment.



#### Advantages:

Pigment clear rate is higher

Pigment removal speed is faster

Refused to stimulate melanin active again risk

Remove melanin colleagues and start repair mechanism

Don't damage the normal tissue



## **SHORT DURATION**

2000ps short duration, painless, no blood.



## **NEW TECHNOLOGY**

Utilize active Q switch picosecond technology.



## STABLE ENERGY

Laser device place on the handle, avoid energy loss.



## **HIGH FREQUENCY**

1-20HZ high frequency, suitable fo more treatments.



Use active Q-Switch Picosecond technology Laser device place on handle, no energy loss



## 1064nm & 532nm

Treatment head can be changed precisely

## TREATMENT EFFECT -Eyebrow Tattootattoo-After **Before Before** After - Age Pigment Chloasma-**Before** After **Before** Blue nevus— Nevus Of Ota -**Before Before** After After

#### Company Profile

#### WHY CHOOSE US

#### Professional OEM, ODM service for Ice laser machine

- 1)In stock 12 hours delivery
- 2) Print any color you want for your machine, make it be your and your client's favorite.
- 3) Print your logo on the machine shell and add it to the system as a welcome interface. Make it exclusive in the world.
- 4) Add any language into the machine system, according to you and your client's requirements.
- 5) Add the Remote Rental System into the machine to do lease business.
- 6) Design an exclusive machine shell for you, form your own brand in the market.
- 7) Design a new interface and system of machine, making it most convenient for you and your clients.
- 8) Develop the newest technology to satisfy you and your client's demand.aser machine hair removal made in US hair removal ice hair removal machine laser diode 3 waves beauty machine



\* TUV ISO 13485:2016 newest standard and more strict for production line inspection

#### \* There is few supplier get Medical CE and ISO13485 from TUV in China market now

- \* KM was founded in 2009, and we employ more than 100 highly valued team members. Every year, our team participates in beauty machine exhibitions to actively learn about the demands of the beauty market from the customers and other industry professionals themselves. This allows us added knowledge to design more competitive products
- \* KM has many medical and production certifications around the world, including (TUV) CE, (TUV) ISO 13485:2016, new type machine shell patent certification, Import and export license, Medical equipment production enterprise license and we have been awarded a High-tech enterprise certificate in China.
- \* KM's enterprise vision is customer-oriented and we are dedicated to providing end customers with high quality and cost efficient beauty equipment. KM continues to work hard for customers in order to become a world famous supplier of beauty machines and medical equipment











sales@wfkmdz.com



wfkmdz.com